

— ABLYBIO, Help Your Research



# CD147 (YD20458) Rabbit mAb

货号: **AYD12518**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse liver,Mouse heart
细胞定位	Cell membrane, Photoreceptor inner segment, Cell projection, cilium, photoreceptor outer segment, Endoplasmic reticulum membrane, Basolateral cell membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Predicted to enable cell-cell adhesion mediator activity; signaling receptor activity; and virus receptor activity. Involved in neural retina development; photoreceptor cell maintenance; and spermatogenesis. Located in several cellular components, including acrosomal membrane; photoreceptor inner segment; and photoreceptor outer segment. Is expressed in several structures, including alimentary system; brain; early conceptus; reproductive system; and sensory organ. Orthologous to human BSG (basigin (Ok blood group)).
基因ID	12215
基因名	Bsg
Swiss	P18572 ( <a href="https://www.uniprot.org/uniprotkb/P18572/entry">https://www.uniprot.org/uniprotkb/P18572/entry</a> )
别名	CD147 (YD20458),CD147 (YD20458) Rabbit mAb,Bsg,Basic immunoglobulin superfamily,HT7 antigen,Membrane glycoprotein gp42

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)