

# Phospho-Rb2 p130 (Ser952) (YD18516) Rabbit mAb

货号: **AYD12520**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	<a href="#">WB</a> <a href="#">IHC-P</a>
推荐浓度	
理论分子量	128kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Neuro-2a
细胞定位	Nucleus
纯化	

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation. Recruits and targets histone methyltransferases KMT5B and KMT5C, leading to epigenetic transcriptional repression. Controls histone H4 'Lys-20' trimethylation. Probably acts as a transcription repressor or by recruiting chromatin-modifying enzymes to promoters. Potent inhibitor of E2F-mediated transcription, associates preferentially with E2F5. Binds to cyclins A and E. Binds to and may be involved in the transforming capacity of the adenovirus E1A protein. May act as a tumor suppressor.
基因ID	5934

基因名	RBL2
Swiss	Q08999
别名	Phospho-Rb2 p130 (Ser952) (YD18516)

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)