

Integrin linked ILK (YD12955) Rabbit mAb

货号: **AYD12535**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,Mouse kidney,293T,NIH/3T3,MCF7,A-549,Mouse heart,Mouse lung
细胞定位	Cell junction, focal adhesion, Cell membrane, Cell projection, lamellipodium, Cytoplasm, myofibril, sarcomere, Nucleus, cytoskeleton, microtubule organizing center, centrosome, cell cortex
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein with a kinase-like domain and four ankyrin-like repeats. The encoded protein associates at the cell membrane with the cytoplasmic domain of beta integrins, where it regulates integrin-mediated signal transduction. Activity of this protein is important in the epithelial to mesenchymal transition, and over-expression of this gene is implicated in tumor growth and metastasis. Alternative splicing results in multiple transcript variants.
基因ID	3611
基因名	ILK

Swiss	Q13418
别名	Integrin linked ILK (YD12955)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn