

— ABLYBIO, Help Your Research



Integrin linked ILK (YD12955) Rabbit mAb

货号: **AYD12535**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,Mouse kidney,293T,NIH/3T3,MCF7,A-549,Mouse heart,Mouse lung
细胞定位	Cell junction, focal adhesion, Cell membrane, Cell projection, lamellipodium, Cytoplasm, myofibril, sarcomere, Nucleus, cytoskeleton, microtubule organizing center, centrosome, cell cortex
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein with a kinase-like domain and four ankyrin-like repeats. The encoded protein associates at the cell membrane with the cytoplasmic domain of beta integrins, where it regulates integrin-mediated signal transduction. Activity of this protein is important in the epithelial to mesenchymal transition, and over-expression of this gene is implicated in tumor growth and metastasis. Alternative splicing results in multiple transcript variants.
基因ID	3611
基因名	ILK
Swiss	Q13418 (https://www.uniprot.org/uniprotkb/Q13418/entry)
别名	Integrin linked ILK (YD12955),Integrin linked ILK (YD12955) Rabbit mAb,ILK,ILK-1,ILK-2,Inactive integrin-linked kinase,p59ILK,ILK1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)