

— ABLYBIO, Help Your Research



Survivin (YD12833) Rabbit mAb

货号: **AYD12540**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P ICC/IF FC
推荐浓度	
理论分子量	16kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293F,Jurkat
细胞定位	Cytoplasm, Nucleus, Chromosome, centromere, cytoskeleton, spindle, kinetochore, Midbody
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the inhibitor of apoptosis (IAP) gene family, which encode negative regulatory proteins that prevent apoptotic cell death. IAP family members usually contain multiple baculovirus IAP repeat (BIR) domains, but this gene encodes proteins with only a single BIR domain. The encoded proteins also lack a C-terminus RING finger domain. Gene expression is high during fetal development and in most tumors, yet low in adult tissues. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
基因ID	332
基因名	BIRC5
Swiss	O15392 (https://www.uniprot.org/uniprotkb/O15392/entry)
别名	Survivin (YD12833),Survivin (YD12833) Rabbit mAb,BIRC5,Apoptosis inhibitor 4,Apoptosis inhibitor survivin,API4,IAP4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)