

— ABLYBIO, Help Your Research



ING1 (YD32927) Rabbit mAb

货号: **AYD12551**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	47kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,HT-29,NIH/3T3,Mouse testis,Mouse brain,Mouse spleen,Mouse lung,Rat brain,Rat liver
细胞定位	Nucleus
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a tumor suppressor protein that can induce cell growth arrest and apoptosis. The encoded protein is a nuclear protein that physically interacts with the tumor suppressor protein TP53 and is a component of the p53 signaling pathway. Reduced expression and rearrangement of this gene have been detected in various cancers. Multiple alternatively spliced transcript variants encoding distinct isoforms have been reported.
基因ID	3621
基因名	ING1
Swiss	Q9UK53 (https://www.uniprot.org/uniprotkb/Q9UK53/entry)
别名	ING1 (YD32927),ING1 (YD32927) Rabbit mAb,ING1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)