

CD14 (YD11979) Rabbit mAb

货号: **AYD12556**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7
细胞定位	Cell membrane, Secreted, Membrane raft, Golgi apparatus
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript variants encoding the same protein.
基因ID	929
基因名	CD14
Swiss	P08571
别名	CD14 (YD11979)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn