

Phospho-RAB10 (Thr73) (YD16675) Rabbit mAb

货号: **AYD12615**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	Cytoplasmic vesicle membrane, Golgi apparatus membrane, Golgi apparatus, trans-Golgi network membrane, Endosome membrane, Recycling endosome membrane, Cytoplasmic vesicle, phagosome membrane, Cytoplasm, cytoskeleton, cilium basal body, Endoplasmic reticulum membrane, perinuclear region, Lyosome
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	RAB10 belongs to the RAS (see HRAS; MIM 190020) superfamily of small GTPases. RAB proteins localize to exocytic and endocytic compartments and regulate intracellular vesicle trafficking (Bao et al., 1998 [PubMed 9918381]).[supplied by OMIM, Mar 2009]
基因ID	10890
基因名	RAB10

Swiss	P61026
别名	Phospho-RAB10 (Thr73) (YD16675)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn