

— ABLYBIO, Help Your Research



# PDGF Receptor alpha (YD15993) Rabbit mAb

货号: **AYD12643**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF FC IP
推荐浓度	
理论分子量	123kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C6,Mouse lung
细胞定位	Cell membrane, Cell projection, cilium, Golgi apparatus
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor or family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides. Studies suggest that this gene plays a role in organ development, wound healing, and tumor progression. Mutations in this gene have been associated with idiopathic hypereosinophilic syndrome, somatic and familial gastrointestinal stromal tumors, and a variety of other cancers.
基因ID	5156
基因名	PDGFRA
Swiss	P16234 ( <a href="https://www.uniprot.org/uniprotkb/P16234/entry">https://www.uniprot.org/uniprotkb/P16234/entry</a> )
别名	PDGF Receptor alpha (YD15993),PDGF Receptor alpha (YD15993) Rabbit mAb,PDGFRA,Alpha platelet-derived growth factor receptor,Alpha-type platelet-derived growth factor receptor,CD140 antigen-like family member A,CD140a antigen,Platelet-derived growth factor alpha receptor

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)