

— ABLYBIO, Help Your Research



# Nectin 2 (YD32503) Rabbit mAb

货号: **AYD12664**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	58kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF-7,SW620,A-549,SKOV3
细胞定位	Cell membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
基因ID	5819
基因名	NECTIN2
Swiss	Q92692 ( <a href="https://www.uniprot.org/uniprotkb/Q92692/entry">https://www.uniprot.org/uniprotkb/Q92692/entry</a> )
别名	Nectin 2 (YD32503),Nectin 2 (YD32503) Rabbit mAb,NECTIN2,Herpes virus entry mediator B,Nectin cell adhesion molecule 2,Poliovirus receptor-related protein 2,HVEB,PRR2,PVRL2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)