

CDC37 (YD31725) Rabbit mAb

货号: **AYD12666**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	44kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,HeLa,HepG2,BT474,HT-1080
细胞定位	Cytoplasm
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is highly similar to Cdc 37, a cell division cycle control protein of <i>Saccharomyces cerevisiae</i> . This protein is a molecular chaperone with specific function in cell signal transduction. It has been shown to form complex with Hsp90 and a variety of protein kinases including CDK4, CDK6, SRC, RAF-1, MOK, as well as eIF2 alpha kinases. It is thought to play a critical role in directing Hsp90 to its target kinases.
基因ID	11140
基因名	CDC37
Swiss	Q16543

别名	CDC37 (YD31725)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn