

— ABLYBIO, Help Your Research



Phospho-Caldesmon/CDM (Ser759) (YD32116) Rabbit mAb

货号: AYD12678

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	93kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431
细胞定位	Cytoplasm, cytoskeleton, myofibril, stress fiber
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a calmodulin- and actin-binding protein that plays an essential role in the regulation of smooth muscle and nonmuscle contraction. The conserved domain of this protein possesses the binding activities to Ca(2+)-calmodulin, actin, tropomyosin, myosin, and phospholipids. This protein is a potent inhibitor of the actin-tropomyosin activated myosin MgATPase, and serves as a mediating factor for Ca(2+)-dependent inhibition of smooth muscle contraction. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms.
基因ID	800
基因名	CALD1
Swiss	Q05682 (https://www.uniprot.org/uniprotkb/Q05682/entry)
别名	Phospho-Caldesmon/CDM (Ser759) (YD32116), Phospho-Caldesmon/CDM (Ser759) (YD32116) Rabbit mAb, CALD1, CAD, CDM

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)