



Nup107 (YD15208) Rabbit mAb

货号: **AYD12683**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	106kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,HeLa,293T,Jurkat,HepG2,Mouse brain,Mouse kidney
细胞定位	Nucleus membrane, Nucleus, nuclear pore complex, Chromosome, centromere, kinetochore
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the nucleoporin family. The protein is localized to the nuclear rim and is an essential component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.
------	---

基因ID	57122
基因名	NUP107
Swiss	P57740 (https://www.uniprot.org/uniprotkb/P57740/entry)
别名	Nup107 (YD15208),Nup107 (YD15208) Rabbit mAb,NUP107,107 kDa nucleoporin,Nucleoporin Nup107

产品验证

 AYD12683-wb1

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)