

— ABLYBIO, Help Your Research



# Phospho-EGFR (Tyr1086) (YD31672) Rabbit mAb

货号: **AYD12700**

## 产品信息

|       |   |
|-------|---|
| 反应    | Human   |
| 宿主    | Rabbit  |
| 克隆性   | Monoclonal  |
| 预测反应  |   |
| 应用    | WB ICC/IF IP  |
| 推荐浓度  |   |
| 理论分子量 | 134kDa  |
| 实测分子量 |   |
| 形式    | Liquid  |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.   |
| 偶联物   | Unconjugated  |
| 阳性对照  | Mouse brain,Rat brain   |
| 细胞定位  | Cell membrane, Endoplasmic reticulum membrane, Golgi apparatus membrane, Nucleus membrane, Endosome, Endosome membrane, Nucleus, Secreted |
| 纯化    | 亲和纯化  |

## 抗原信息

|      |  |
|------|--|
| 抗原信息 |  |
|------|--|

## 靶点信息

|       |   |
|-------|---|
| 研究背景  | The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer. |
| 基因ID  | 1956  |
| 基因名   | EGFR  |
| Swiss | P00533 ( <a href="https://www.uniprot.org/uniprotkb/P00533/entry">https://www.uniprot.org/uniprotkb/P00533/entry</a> )  |
| 别名    | Phospho-EGFR (Tyr1086) (YD31672), Phospho-EGFR (Tyr1086) (YD31672) Rabbit mAb, EGFR, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1, ERBB, ERBB1, HER1  |

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))