

— ABLYBIO, Help Your Research



# DAG1 (YD33473) Rabbit mAb

货号: **AYD12706**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	97kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U2OS,SW480,HeLa,mouse skeletal muscle,mouse lung,mouse heart,rat liver,rat lung
细胞定位	Secreted, extracellular space, extracellular matrix, basement membrane, Synapse, Cell membrane, Cyto plasm, cytoskeleton, Nucleus, nucleoplasm, sarcolemma, Postsynaptic cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes dystroglycan, a central component of dystrophin-glycoprotein complex that links the extracellular matrix and the cytoskeleton in the skeletal muscle. The encoded preproprotein undergoes O- and N-glycosylation, and proteolytic processing to generate alpha and beta subunits. Certain mutations in this gene are known to cause distinct forms of muscular dystrophy. Alternative splicing results in multiple transcript variants, all encoding the same protein.
基因ID	1605
基因名	DAG1
Swiss	Q14118 ( <a href="https://www.uniprot.org/uniprotkb/Q14118/entry">https://www.uniprot.org/uniprotkb/Q14118/entry</a> )
别名	DAG1 (YD33473), DAG1 (YD33473) Rabbit mAb, DAG1, Dystroglycan, Dystrophin-associated glycoprotein 1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)