

— ABLYBIO, Help Your Research



# Phospho-LRRK2 (Ser955) (YD31380) Rabbit mAb

货号: AYD12712

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	286kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cytoplasmic vesicle, Perikaryon, Golgi apparatus membrane, Cell projection, axon, dendrite, Endoplasmic reticulum membrane, secretory vesicle, synaptic vesicle membrane, Endosome, Lysosome, Mitochondrion outer membrane, Cytoplasm, cytoskeleton, phagosome
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene is a member of the leucine-rich repeat kinase family and encodes a protein with an ankryin repeat region, a leucine-rich repeat (LRR) domain, a kinase domain, a DFG-like motif, a RAS domain, a GTPase domain, a MLK-like domain, and a WD40 domain. The protein is present largely in the cytoplasm but also associates with the mitochondrial outer membrane. Mutations in this gene have been associated with Parkinson disease-8.
基因ID	120892
基因名	LRRK2
Swiss	Q5S007 ( <a href="https://www.uniprot.org/uniprotkb/Q5S007/entry">https://www.uniprot.org/uniprotkb/Q5S007/entry</a> )
别名	Phospho-LRRK2 (Ser955) (YD31380),Phospho-LRRK2 (Ser955) (YD31380) Rabbit mAb,LRRK2,Dardarin,PA RK8

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)