

— ABLYBIO, Help Your Research



# Haptoglobin (YD20954) Rabbit mAb

货号: **AYD12749**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,MCF7,HL-60,Mouse liver,Mouse lung,Mouse spleen,Rat lung,Rat spleen
细胞定位	Secreted
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a preproprotein, which is processed to yield both alpha and beta chains, which subsequently combine as a tetramer to produce haptoglobin. Haptoglobin functions to bind free plasma hemoglobin, which allows degradative enzymes to gain access to the hemoglobin, while at the same time preventing loss of iron through the kidneys and protecting the kidneys from damage by hemoglobin. Mutations in this gene and/or its regulatory regions cause ahaptoglobinemia or hypohaptoglobinemia. This gene has also been linked to diabetic nephropathy, the incidence of coronary artery disease in type 1 diabetes, Crohn's disease, inflammatory disease behavior, primary sclerosing cholangitis, susceptibility to idiopathic Parkinson's disease, and a reduced incidence of Plasmodium falciparum malaria. The protein encoded also exhibits antimicrobial activity against bacteria. A similar duplicated gene is located next to this gene on chromosome 16. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3240
基因名	HP
Swiss	P00738 ( <a href="https://www.uniprot.org/uniprotkb/P00738/entry">https://www.uniprot.org/uniprotkb/P00738/entry</a> )
别名	Haptoglobin (YD20954),Haptoglobin (YD20954) Rabbit mAb,HP,Zonulin

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)