

— ABLYBIO, Help Your Research



# CD47 (YD20943) Rabbit mAb

货号: **AYD12754**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	ICC/IF FC
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse brain,Rat brain
细胞定位	Cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may play a role in membrane transport and signal transduction. This gene has broad tissue distribution, and is reduced in expression on Rh erythrocytes. Alternatively spliced transcript variants have been found for this gene.
基因ID	961
基因名	CD47
Swiss	Q08722 ( <a href="https://www.uniprot.org/uniprotkb/Q08722/entry">https://www.uniprot.org/uniprotkb/Q08722/entry</a> )
别名	CD47 (YD20943),CD47 (YD20943) Rabbit mAb,CD47,Antigenic surface determinant protein OA3,Integrin-associated protein,Protein MER6,MER6

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)