

ABCB11/BSEP (YD20810) Rabbit mAb

货号: **AYD12767**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P ICC/IF FC
推荐浓度	
理论分子量	146kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Apical cell membrane, Recycling endosome membrane, Endosome, Cell membrane
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN2, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by this gene is the major canalicular bile salt export pump in man. Mutations in this gene cause a form of progressive familial intrahepatic cholestases which are a group of inherited disorders with severe cholestatic liver disease from early infancy.
基因ID	8647
基因名	ABCB11

Swiss	O95342
别名	ABCB11/BSEP (YD20810)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn