

MAP2 (YD20476) Rabbit mAb

货号: **AYD12803**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	200kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,PC-12
细胞定位	Cytoplasm, cytoskeleton, Cell projection, dendrite
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described.
基因ID	4133
基因名	MAP2
Swiss	P11137

别名	MAP2 (YD20476)
----	----------------

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn