

— ABLYBIO, Help Your Research



KDM5A/Jarid1A/RBBP2 (YD20451) Rabbit mAb

货号: **AYD12806**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	192kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,MCF7
细胞定位	Nucleus, nucleolus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the Jumonji, AT-rich interactive domain 1 (JARID1) histone demethylase protein family. The encoded protein plays a role in gene regulation through the histone code by specifically demethylating lysine 4 of histone H3. The encoded protein interacts with many other proteins, including retinoblastoma protein, and is implicated in the transcriptional regulation of Hox genes and cytokines. This gene may play a role in tumor progression.
基因ID	5927
基因名	KDM5A
Swiss	P29375 (https://www.uniprot.org/uniprotkb/P29375/entry)
别名	KDM5A/Jarid1A/RBBP2 (YD20451),KDM5A/Jarid1A/RBBP2 (YD20451) Rabbit mAb,KDM5A,Histone demethylase JARID1A,Jumonji/ARID domain-containing protein 1A,Retinoblastoma-binding protein 2,[histone H3]-trimethyl-L-lysine(4) demethylase 5A,JARID1A,RBBP2,RBP2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)