

— ABLYBIO, Help Your Research



# hHR23A (YD33399) Rabbit mAb

货号: **AYD12834**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is one of two human homologs of <i>Saccharomyces cerevisiae</i> Rad23, a protein involved in nucleotide excision repair. Proteins in this family have a modular domain structure consisting of an ubiquitin-like domain (Ubl), ubiquitin-associated domain 1 (Uba1), XPC-binding domain and Uba2. The protein encoded by this gene plays an important role in nucleotide excision repair and also in delivery of polyubiquitinated proteins to the proteasome. Alternative splicing results in multiple transcript variants encoding multiple isoforms.
基因ID	5886
基因名	RAD23A
Swiss	P54725 ( <a href="https://www.uniprot.org/uniprotkb/P54725/entry">https://www.uniprot.org/uniprotkb/P54725/entry</a> )
别名	hHR23A (YD33399),hHR23A (YD33399) Rabbit mAb,RAD23A,UV excision repair protein RAD23 homolog A

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)