

— ABLYBIO, Help Your Research



# CD63 (YD19296) Rabbit mAb

货号: **AYD12861**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P ICC/IF FC
推荐浓度	
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375
细胞定位	Cell membrane, Lysosome membrane, Late endosome membrane, Endosome, multivesicular body, Melanosome, Secreted, extracellular exosome, Cell surface
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The encoded protein is a cell surface glycoprotein that is known to complex with integrins. It may function as a blood platelet activation marker. Deficiency of this protein is associated with Hermansky-Pudlak syndrome. Also this gene has been associated with tumor progression. Alternative splicing results in multiple transcript variants encoding different protein isoforms.
基因ID	967
基因名	CD63
Swiss	P08962 ( <a href="https://www.uniprot.org/uniprotkb/P08962/entry">https://www.uniprot.org/uniprotkb/P08962/entry</a> )
别名	CD63 (YD19296),CD63 (YD19296) Rabbit mAb,CD63,Granulophysin,Lysosomal-associated membrane protein 3,Lysosome integral membrane protein 1,Melanoma-associated antigen ME491,OMA81H,Ocular melanoma-associated antigen,Tetraspanin-30,MLA1,TSPAN30

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)