

— ABLYBIO, Help Your Research



# LAMP2 (YD20923) Rabbit mAb

货号: **AYD12870**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,HepG2,Mouse liver
细胞定位	Lysosome membrane, Endosome membrane, Cell membrane, Cytoplasmic vesicle, autophagosome membrane
纯化	亲和纯化




## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins.
基因ID	3920
基因名	LAMP2
Swiss	P13473 ( <a href="https://www.uniprot.org/uniprotkb/P13473/entry">https://www.uniprot.org/uniprotkb/P13473/entry</a> )
别名	LAMP2 (YD20923),LAMP2 (YD20923) Rabbit mAb,LAMP2,CD107 antigen-like family member B,LGP-96

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)