

CD83 (YD18735) Rabbit mAb

货号: **AYD12877**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse intestine
细胞定位	Membrane
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene.
基因ID	9308
基因名	CD83
Swiss	Q01151

别名	CD83 (YD18735)
----	----------------

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn