

— ABLYBIO, Help Your Research



Flt3/CD135 (YD18608) Rabbit mAb

货号: **AYD12885**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	FC
推荐浓度	
理论分子量	113kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Membrane, Endoplasmic reticulum lumen
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a class III receptor tyrosine kinase that regulates hematopoiesis. This receptor is activated by binding of the fms-related tyrosine kinase 3 ligand to the extracellular domain, which induces homodimer formation in the plasma membrane leading to autophosphorylation of the receptor. The activated receptor kinase subsequently phosphorylates and activates multiple cytoplasmic effector molecules in pathways involved in apoptosis, proliferation, and differentiation of hematopoietic cells in bone marrow. Mutations that result in the constitutive activation of this receptor result in acute myeloid leukemia and acute lymphoblastic leukemia.
基因ID	2322
基因名	FLT3
Swiss	P36888 (https://www.uniprot.org/uniprotkb/P36888/entry)
别名	Flt3/CD135 (YD18608), Flt3/CD135 (YD18608) Rabbit mAb, FLT3, FL cytokine receptor, Fetal liver kinase-2, Fms-like tyrosine kinase 3, Stem cell tyrosine kinase 1, CD135, FLK2, STK1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)