

— ABLYBIO, Help Your Research



# Cytokeratin 13 (YD34223) Rabbit mAb

货号: **AYD12952**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse intestine,Rat thymus
细胞定位	
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the supra basal layers of non-cornified stratified epithelia. Mutations in this gene and keratin 4 have been associated with the autosomal dominant disorder White Sponge Nevus. The type I cytokeratins are clustered in a region of chromosome 17q21.2. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been described.
基因ID	3860
基因名	KRT13
Swiss	P13646 ( <a href="https://www.uniprot.org/uniprotkb/P13646/entry">https://www.uniprot.org/uniprotkb/P13646/entry</a> )
别名	Cytokeratin 13 (YD34223),Cytokeratin 13 (YD34223) Rabbit mAb,KRT13,Cytokeratin-13,Keratin-13

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)