

— ABLYBIO, Help Your Research



Desmoglein 3/PVA (YD34220) Rabbit mAb

货号: **AYD12953**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	108kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	A431
细胞定位	Cell membrane, Cell junction, desmosome, Cytoplasm, tight junction
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the desmoglein family and cadherin cell adhesion molecule superfamily of proteins. Desmogleins are calcium-binding transmembrane glycoprotein components of desmosomes, cell-cell junctions between epithelial, myocardial, and other cell types. The encoded preproprotein is proteolytically processed to generate the mature glycoprotein. This gene is present in a gene cluster with other desmoglein gene family members on chromosome 18. The encoded protein has been identified as the autoantigen of the autoimmune blistering disease pemphigus vulgaris.
基因ID	1830
基因名	DSG3
Swiss	P32926 (https://www.uniprot.org/uniprotkb/P32926/entry)
别名	Desmoglein 3/PVA (YD34220),Desmoglein 3/PVA (YD34220) Rabbit mAb,DSG3,130 kDa pemphigus vulgaris antigen,Cadherin family member 6,CDHF6

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)