

— ABLYBIO, Help Your Research



# Ube2N/Ubc13 (YD34009) Rabbit mAb

货号: **AYD13001**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,B cells,Mouse brain,Mouse spleen,Rat brain,Rat testis
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Studies in mouse suggest that this protein plays a role in DNA postreplication repair.
基因ID	7334
基因名	UBE2N
Swiss	P61088 ( <a href="https://www.uniprot.org/uniprotkb/P61088/entry">https://www.uniprot.org/uniprotkb/P61088/entry</a> )
别名	Ube2N/Ubc13 (YD34009), Ube2N/Ubc13 (YD34009) Rabbit mAb, UBE2N, Bendless-like ubiquitin-conjugating enzyme, E2 ubiquitin-conjugating enzyme N, Ubc13, UbcH13, Ubiquitin carrier protein N, Ubiquitin-protein ligase N, BLU

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)