

— ABLYBIO, Help Your Research



# ATP synthase C (YD11540) Rabbit mAb

货号: **AYD13011**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	14kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,Mouse heart,Mouse small intestine,Mouse liver,Rat heart,Rat brain
细胞定位	Mitochondrion membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel seems to have nine subunits (a, b, c, d, e, f, g, F6 and 8). This gene is one of three genes that encode subunit c of the proton channel. Each of the three genes have distinct mitochondrial import sequences but encode the identical mature protein. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]
基因ID	516
基因名	ATP5MC1
Swiss	P05496 ( <a href="https://www.uniprot.org/uniprotkb/P05496/entry">https://www.uniprot.org/uniprotkb/P05496/entry</a> )
别名	ATP synthase C (YD11540),ATP synthase C (YD11540) Rabbit mAb,ATP5MC1,ATP synthase lipid-binding protein,ATP synthase membrane subunit c locus 1,ATP synthase proteolipid P1,ATP synthase proton-transporting mitochondrial F(0) complex subunit C1,ATPase protein 9,ATPase subunit c

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)