

— ABLYBIO, Help Your Research



# Carbonic Anhydrase 1/CA1 (YD33770) Rabbit mAb

货号: AYD13046

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse liver,Mouse heart,Mouse brain
细胞定位	Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. This CA1 gene is closely linked to the CA2 and CA3 genes on chromosome 8. It encodes a cytosolic protein that is found at the highest level in erythrocytes. Allelic variants of this gene have been described in some populations. Alternative splicing and the use of alternative promoters results in multiple transcript variants.
基因ID	759
基因名	CA1
Swiss	P00915 ( <a href="https://www.uniprot.org/uniprotkb/P00915/entry">https://www.uniprot.org/uniprotkb/P00915/entry</a> )
别名	Carbonic Anhydrase 1/CA1 (YD33770),Carbonic Anhydrase 1/CA1 (YD33770) Rabbit mAb,CA1,Carbonate dehydratase I,Carbonic anhydrase B,Carbonic anhydrase I,Cyanamide hydratase CA1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)