

— ABLYBIO, Help Your Research



Ficolin 2 (YD33758) Rabbit mAb

货号: **AYD13050**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	34kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mosue spleen,mouse lung,rat spleen
细胞定位	Secreted, Cell surface
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The product of this gene belongs to the ficolin family of proteins. This family is characterized by the presence of a leader peptide, a short N-terminal segment, followed by a collagen-like region, and a C-terminal fibrinogen-like domain. This gene is predominantly expressed in the liver, and has been shown to have carbohydrate binding and opsonic activities. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	2220
基因名	FCN2
Swiss	Q15485 (https://www.uniprot.org/uniprotkb/Q15485/entry)
别名	Ficolin 2 (YD33758),Ficolin 2 (YD33758) Rabbit mAb,FCN2,37 kDa elastin-binding protein,Collagen/fibrinogen domain-containing protein 2,Ficolin-B,Ficolin-beta,Hucolin,L-ficolin,Serum lectin p35,FCNL

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)