

— ABLYBIO, Help Your Research



# Carbonic Anhydrase 9/CA9 (YD33749) Rabbit mAb

货号: **AYD13052**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Nucleus, nucleolus, Cell membrane, Cell projection, microvillus membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. CA IX is a transmembrane protein and is one of only two tumor-associated carbonic anhydrase isoenzymes known. It is expressed in all clear-cell renal cell carcinoma, but is not detected in normal kidney or most other normal tissues. It may be involved in cell proliferation and transformation. This gene was mapped to 17q21.2 by fluorescence in situ hybridization, however, radiation hybrid mapping localized it to 9p13-p12.
基因ID	768
基因名	CA9
Swiss	Q16790 ( <a href="https://www.uniprot.org/uniprotkb/Q16790/entry">https://www.uniprot.org/uniprotkb/Q16790/entry</a> )
别名	Carbonic Anhydrase 9/CA9 (YD33749),Carbonic Anhydrase 9/CA9 (YD33749) Rabbit mAb,CA9,Carbonate dehydratase IX,Carbonic anhydrase IX,Membrane antigen MN,P54/58N,Renal cell carcinoma-associated antigen G250,pMW1,G250,MN

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)