

eIF4B (YD33486) Rabbit mAb

货号: **AYD13058**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat skeletal muscle
细胞定位	
纯化	




抗原信息

抗原信息	
------	--

靶点信息

研究背景	Required for the binding of mRNA to ribosomes. Functions in close association with EIF4-F and EIF4-A. Binds near the 5'-terminal cap of mRNA in presence of EIF-4F and ATP. Promotes the ATPase activity and the ATP-dependent RNA unwinding activity of both EIF4-A and EIF4-F.
基因ID	1975
基因名	EIF4B
Swiss	P23588
别名	eIF4B (YD33486)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn