

VCP (YD33642) Rabbit mAb

货号: **AYD13081**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	89kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,SH-SY5Y,Mouse brain,Mouse spleen,Mouse kidney,Rat liver,Rat brain
细胞定位	Cytoplasm, cytosol, Endoplasmic reticulum, Nucleus, Stress granule
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of a family that includes putative ATP-binding proteins involved in vesicle transport and fusion, 26S proteasome function, and assembly of peroxisomes. This protein, as a structural protein, is associated with clathrin, and heat-shock protein Hsc70, to form a complex. It has been implicated in a number of cellular events that are regulated during mitosis, including homotypic membrane fusion, spindle pole body function, and ubiquitin-dependent protein degradation.
基因ID	7415
基因名	VCP
Swiss	P55072

别名	VCP (YD33642)
----	---------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn