

— ABLYBIO, Help Your Research



FES (YD33640) Rabbit mAb

货号: **AYD13082**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	93kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse brain
细胞定位	Cytoplasm, cytosol, cytoskeleton, Cell membrane, Cytoplasmic vesicle, Golgi apparatus, Cell junction, focal adhesion
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes the human cellular counterpart of a feline sarcoma retrovirus protein with transforming capabilities. The gene product has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular transformation. Its chromosomal location has linked it to a specific translocation event identified in patients with acute promyelocytic leukemia but it is also involved in normal hematopoiesis as well as growth factor and cytokine receptor signaling. Alternative splicing results in multiple variants encoding different isoforms.
基因ID	2242
基因名	FES
Swiss	P07332 (https://www.uniprot.org/uniprotkb/P07332/entry)
别名	FES (YD33640),FES (YD33640) Rabbit mAb,FES,Feline sarcoma/Fujinami avian sarcoma oncogene homolog,Proto-oncogene c-Fes,Proto-oncogene c-Fps,p93c-fes,FPS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)