

— ABLYBIO, Help Your Research



PACAP (YD33600) Rabbit mAb

货号: **AYD13090**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	Endoplasmic reticulum lumen, Secreted, Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Associates with immunoglobulin M (IgM heavy and light chains and promotes IgM assembly and secretion. May exert its effect by acting as a molecular chaperone or as an oxidoreductase as it displays a low level of oxidoreductase activity (By similarity. Isoform 2 may be involved in regulation of apoptosis. Helps to diversify peripheral B-cell functions by regulating Ca ²⁺ stores, antibody secretion and integrin activation.
基因ID	51237
基因名	MZB1
Swiss	Q8WU39 (https://www.uniprot.org/uniprotkb/Q8WU39/entry)
别名	PACAP (YD33600),PACAP (YD33600) Rabbit mAb,MZB1,Mesenteric estrogen-dependent adipose 7,Plasma cell-induced resident endoplasmic reticulum protein,Proapoptotic caspase adapter protein,MEDA7,PACAP

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)