

UCH37 (YD33598) Rabbit mAb

货号: **AYD13091**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	38kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Neuro-2a,PC-12
细胞定位	Cytoplasm, Nucleus
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	Uch37 is one of the deubiquitinating enzymes of the human 26S proteasome, belongs to the Ub C-terminal hydrolases. Uch37 was found to bind both Adrm1 and S5a subunits of the 26S proteasome. Uch37 likely contributes to the Ub chain trimming activity of the 26S proteasome that shortens a polyUb chain from the distal end of the chain.
基因ID	51377
基因名	UCHL5
Swiss	Q9Y5K5

别名	UCH37 (YD33598)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn