

— ABLYBIO, Help Your Research



Argonaute-2 (YD12085) Rabbit mAb

货号: **AYD13092**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	97kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,293T,HeLa
细胞定位	Cytoplasm, P-body, Nucleus
纯化	亲和纯化

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, and contains a PAZ domain and a PIWI domain. It may interact with di cer1 and play a role in short-interfering-RNA-mediated gene silencing. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	27161
基因名	AGO2
Swiss	Q9UKV8 (https://www.uniprot.org/uniprotkb/Q9UKV8/entry)
别名	Argonaute-2 (YD12085),Argonaute-2 (YD12085) Rabbit mAb,AGO2,Argonaute RISC catalytic component 2,Eukaryotic translation initiation factor 2C 2,PAZ Piwi domain protein,Protein slicer,EIF2C2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)