

MPP1 (YD33403) Rabbit mAb

货号: **AYD13110**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane, Cell projection, stereocilium
纯化	



抗原信息

抗原信息	
------	--

靶点信息

研究背景	Essential regulator of neutrophil polarity. Regulates neutrophil polarization by regulating AKT1 phosphorylation through a mechanism that is independent of PIK3CG activity (By similarity)
基因ID	4354
基因名	MPP1
Swiss	Q00013
别名	MPP1 (YD33403)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn