

— ABLYBIO, Help Your Research



# Vitronectin (YD18920) Rabbit mAb

货号: **AYD13111**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human serum,Rat serum
细胞定位	Secreted, extracellular space, Parasitophorous vacuole
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.
------	--

基因ID	7448
基因名	VTN
Swiss	P04004 ( <a href="https://www.uniprot.org/uniprotkb/P04004/entry">https://www.uniprot.org/uniprotkb/P04004/entry</a> )
别名	Vitronectin (YD18920), Vitronectin (YD18920) Rabbit mAb, VTN, S-protein, Serum-spreading factor, V75

## 产品验证

Human plasma

kDa

180 -

130 -

100 -

70 -

55 -

40 -

35 -

25 -

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)