

Vitronectin (YD18920) Rabbit mAb

货号: **AYD13111**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human serum,Rat serum
细胞定位	Secreted, extracellular space, Parasitophorous vacuole
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.
基因ID	7448
基因名	VTN
Swiss	P04004

别名	Vitronectin (YD18920)
----	-----------------------

产品验证

Human plasma

kDa
180 -
130 -
100 -
70 -
55 -
40 -
35 -
25 -



实验步骤

访问官网浏览详情: www.ablybio.cn