

— ABLYBIO, Help Your Research



IL-36 alpha/IL-1F6 (YD33539) Rabbit mAb

货号: **AYD13116**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,BT-474,SW480,A-431,HeLa,Mouse skin,Mouse heart,Rat liver
细胞定位	Cytoplasm, Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a cytokine that can activate NF-kappa-B and MAPK signaling pathways to generate an inflammatory response. The encoded protein functions primarily in skin and demonstrates increased expression in psoriasis. In addition, decreased expression of this gene has been linked to a poor prognosis in both hepatocellular carcinoma and colorectal cancer patients.
------	---

基因ID	27179
基因名	IL36A
Swiss	Q9UHA7 (https://www.uniprot.org/uniprotkb/Q9UHA7/entry)
别名	IL-36 alpha/IL-1F6 (YD33539),IL-36 alpha/IL-1F6 (YD33539) Rabbit mAb,IL36A,FIL1 epsilon,Interleukin-1 epsilon,Interleukin-1 family member 6,FIL1E,IL1E,IL1F6

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)