

— ABLYBIO, Help Your Research



# Myelin oligodendrocyte glycoprotein (YD33414) Rabbit mAb

货号: AYD13154

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The product of this gene is a membrane protein expressed on the oligodendrocyte cell surface and the outermost surface of myelin sheaths. Due to this localization, it is a primary target antigen involved in immune-mediated demyelination. This protein may be involved in completion and maintenance of the myelin sheath and in cell-cell communication. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	4340
基因名	MOG
Swiss	Q16653 ( <a href="https://www.uniprot.org/uniprotkb/Q16653/entry">https://www.uniprot.org/uniprotkb/Q16653/entry</a> )
别名	Myelin oligodendrocyte glycoprotein (YD33414), Myelin oligodendrocyte glycoprotein (YD33414) Rabbit mAb, MOG

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)