

— ABLYBIO, Help Your Research



# Prothrombin (YD14409) Rabbit mAb

货号: **AYD13156**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Thrombin, which cleaves bonds after Arg and Lys, converts fibrinogen to fibrin and activates factors V, VI, VIII, XIII, and, in complex with thrombomodulin, protein C. Functions in blood homeostasis, inflammation and wound healing (By similarity). Activates coagulation factor XI (F11); activation is promoted by the contact with negatively charged surfaces (By similarity). Triggers the production of pro-inflammatory cytokines, such as MCP-1/CCL2 and IL8/CXCL8, in endothelial cells (By similarity)
基因ID	3134
基因名	F2
Swiss	P19221 ( <a href="https://www.uniprot.org/uniprotkb/P19221/entry">https://www.uniprot.org/uniprotkb/P19221/entry</a> )
别名	Prothrombin (YD14409),Prothrombin (YD14409) Rabbit mAb,F2,Coagulation factor II,Cf2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)