

— ABLYBIO, Help Your Research



Prothrombin (YD33407) Rabbit mAb

货号: **AYD13158**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Secreted, extracellular space
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Coagulation factor II is proteolytically cleaved to form thrombin in the first step of the coagulation cascade which ultimately results in the stemming of blood loss. F2 also plays a role in maintaining vascular integrity during development and postnatal life. Peptides derived from the C-terminus of this protein have an antimicrobial activity against E. coli and P. aeruginosa. Mutations in F2 lead to various forms of thrombosis and dysprothrombinemia. Alternative splicing results in multiple transcript variants.
基因ID	2147
基因名	F2
Swiss	P00734 (https://www.uniprot.org/uniprotkb/P00734/entry)
别名	Prothrombin (YD33407),Prothrombin (YD33407) Rabbit mAb,F2,Coagulation factor II

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)