

— ABLYBIO, Help Your Research



# TCP1 alpha/CCTA (YD33291) Rabbit mAb

货号: **AYD13207**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF IP
推荐浓度	
理论分子量	60kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,PC-3,HeLa,Jurkat
细胞定位	Cytoplasm, cytosol, cytoskeleton, microtubule organizing center, centrosome
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a molecular chaperone that is a member of the chaperonin containing TCP1 complex (CCT), also known as the TCP1 ring complex (TRiC). This complex consists of two identical stacked rings, each containing eight different proteins. Unfolded polypeptides enter the central cavity of the complex and are folded in an ATP-dependent manner. The complex folds various proteins, including actin and tubulin. Alternate transcriptional splice variants of this gene, encoding different isoforms, have been characterized. In addition, three pseudogenes that appear to be derived from this gene have been found.
基因ID	6950
基因名	TCP1
Swiss	P17987 ( <a href="https://www.uniprot.org/uniprotkb/P17987/entry">https://www.uniprot.org/uniprotkb/P17987/entry</a> )
别名	TCP1 alpha/CCTA (YD33291), TCP1 alpha/CCTA (YD33291) Rabbit mAb, TCP1, CCT-alpha, Chaperonin containing T-complex polypeptide 1 subunit 1, CCT1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)