

— ABLYBIO, Help Your Research



APOA4/Apo-AIV (YD33289) Rabbit mAb

货号: **AYD13208**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse small intestine,Rat liver
细胞定位	Secreted
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Apolipoprotein (apo) A-IV gene contains 3 exons separated by two introns. A sequence polymorphism has been identified in the 3'UTR of the third exon. The primary translation product is a 396-residue preprotein which after proteolytic processing is secreted its primary site of synthesis, the intestine, in association with chylomicron particles. Although its precise function is not known, apo A-IV is a potent activator of lecithin-cholesterol acyltransferase in vitro.
基因ID	337
基因名	APOA4
Swiss	P06727 (https://www.uniprot.org/uniprotkb/P06727/entry)
别名	APOA4/Apo-AIV (YD33289),APOA4/Apo-AIV (YD33289) Rabbit mAb,APOA4,Apolipoprotein A4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)