

— ABLYBIO, Help Your Research



FKBP6 (YD33279) Rabbit mAb

货号: **AYD13218**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,Mouse heart,Mouse liver,Rat testis,Rat heart
细胞定位	Cytoplasm, Nucleus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a cis-trans peptidyl-prolyl isomerase that may function in immunoregulation and basic cellular processes involving protein folding and trafficking. This gene is located in a chromosomal region that is deleted in Williams-Beuren syndrome. Defects in this gene may cause male infertility. There are multiple pseudogenes for this gene located nearby on chromosome 7. Alternative splicing results in multiple transcript variants.
基因ID	8468
基因名	FKBP6
Swiss	O75344 (https://www.uniprot.org/uniprotkb/O75344/entry)
别名	FKBP6 (YD33279),FKBP6 (YD33279) Rabbit mAb,FKBP6,36 kDa FK506-binding protein,FK506-binding protein 6,Immunophilin FKBP36,FKBP36

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)