

— ABLYBIO, Help Your Research



# Cytokeratin 7 (YD11545) Rabbit mAb

货号: **AYD13225**

## 产品信息

|       |   |
|-------|---|
| 反应    | Human   |
| 宿主    | Rabbit  |
| 克隆性   | Monoclonal  |
| 预测反应  |   |
| 应用    | WB IP   |
| 推荐浓度  |   |
| 理论分子量 | 51kDa   |
| 实测分子量 |   |
| 形式    | Liquid  |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated  |
| 阳性对照  | HeLa  |
| 细胞定位  | Cytoplasm   |
| 纯化    | 亲和纯化  |


## 抗原信息

|      |  |
|------|--|
| 抗原信息 |  |
|------|--|

## 靶点信息

|       |  |
|-------|--|
| 研究背景  | The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the simple epithelia lining the cavities of the internal organs and in the gland ducts and blood vessels. The genes encoding the type II cytokeratins are clustered in a region of chromosome 12q12-q13. Alternative splicing may result in several transcript variants; however, not all variants have been fully described. |
| 基因ID  | 3855   |
| 基因名   | KRT7   |
| Swiss | P08729 ( <a href="https://www.uniprot.org/uniprotkb/P08729/entry">https://www.uniprot.org/uniprotkb/P08729/entry</a> )   |
| 别名    | Cytokeratin 7 (YD11545),Cytokeratin 7 (YD11545) Rabbit mAb,KRT7,Cytokeratin-7,Keratin-7,Sarcolectin,Type-II keratin Kb7,SCL  |

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)